

Current Status: Active PolicyStat ID: 7428064



 Origination:
 04/1999

 Effective:
 01/2020

 Last Approved:
 01/2020

 Last Revised:
 01/2020

 Next Review:
 01/2023

Owner: Roshelle Satterthwait: Dir

Perioperative Svcs

Policy Area: Sterile Processing

References:

Applicability: WA - Kadlec Regional Medical

Center

# Instrument Care and Cleaning Through Decontamination, 32.11.02

**Document Type:Policy** 

### **DESCRIPTION:**

This document provides guidelines to assist Sterile Processing staff in decontaminating and preparing surgical instruments and powered equipment for terminal sterilization and disinfection.

## **PURPOSE:**

To ensure instruments, equipment and case carts are efficiently and effectively decontaminated.

## SCOPE:

Begins with the receipt of the contaminated item ending with the start of the wash cycle.

#### **GENERAL:**

- a. All items received into Decontamination shall be considered contaminated and require cleaning.
- b. All items transported to the SPD Decontamination area shall be covered and contained to prevent exposure of patients or personnel to blood borne pathogens and other potentially infectious organisms.
- c. Upon entrance to Decontamination, SPD Techs mush wash their hands.
- d. SPD Staff handling contaminated instruments and equipment shall wear appropriate personal protective equipment (PPE) including face mask, hair cover, boot covers, an impervious gown, and gloves.
- e. SPD Staff prepares decontamination room before and after their shift.
  - 1. Disinfects counters, sinks, and all surface areas.
  - 2. Dispose of any trash.
  - 3. Restock any supplies that may be needed.
  - 4. Clean screens of washers / cart wash
  - 5. Fill and degas ultrasonics
- f. Items received shall be decontaminated according to manufacturers' instructions (IFU) regarding types of

cleaning agents used (Enzymatic preparations, Detergents) and types of cleaning methods (Manual, Automated).

- g. Instruments shall be arranged into trays for optimal cleaning and drying
- h. Washer racks will be loaded to facilitate the most surface area contact possible.
- Special precautions shall be taken for reprocessing ophthalmic surgical instruments to prevent contaminants being introduced into the eye during surgery (refer to Standards of Work for Decontamination for more information)

(For further information on Decontamination processes refer to the Standards of Work)

# **Attachments**

No Attachments

# **Approval Signatures**

Approver	Date
Kirk Harper: CNO	01/2020
Heather Shipman: Executive Assistant	01/2020
Roshelle Satterthwait: Dir Perioperative Svcs	01/2020

# **Applicability**

WA - Kadlec Regional Medical Center