



Send Out Collection Materials



Quest Viral Cervical Transport Media



Quest Nasopharyngeal Viral Transport Media


Pink and Gray Para-paks



Black Top Total-Fix



Gene Expert Materials




FRONT VIEW



BACK VIEW

**CT-NG Gene Xpert
Material for Urine**
For CT/NG, and Trichomonas PCR
testing.




FRONT VIEW



BACK VIEW

**CT-NG, Trichomonas, and
MVP (Vaginal Panel) GeneXpert
Material for
Female (Swab)**

Viral and Fungal Materials



**Viral Transport Media
For: COVID, Influenza, RSV,
Respiratory Panel PCR, and
Sendout Viral Culture testing.**



**Orange Para-Pak for
Stool PCR Testing.**

Swabs


Blue Swab

For Aerobic and Anaerobic cultures.



ESwab

For Aerobic and Anaerobic cultures and Stool PCR Testing.



Red Swab

For: Strep A testing, Aerobic cultures only, and MRSA PCR Nasal tests.



BD Flocked Swab

Influenza, RSV, & COVID (at CCH only)

For COVID swab is placed in a transport.



For COVID only

